

Powerful!
Mobile X-Ray System

AceMobil-510

Powerful 100 mA output with Small Build



*Arm Free Stop function at any position
High Frequency Inverter system
Maximum mA Current 100 mA
Wireless & Wire Remote Exposure
Light Weight & Low Centroid Design*





Powerful! Mobile X-Ray System

AceMobil-510

Main Item	Sub Item	Descriptions
Generator Type		High Frequency Inverter Control System
Output	Adjustable Tube Voltage	40 kV~110 kV (1 kV Step)
	mAs Range	1.0 mAs ~ 200 mAs
	Max output power	5.1 kW @ 60 kV 85 mA at Input 220 V AC
X-Ray tube	Focal Spot	1.8 mm
	Max. Tube Voltage/Current	More Than 110 kV / 100 mA
	Anode Type	Stationary Anode
	Anode Heat Storage Capacity	More than 30 kJ
Collimator	Target Angle	15°
	Ramp Illumination	More than 160 lux (24 V / 150 W)
	SID Scale	2 m
AC input	Voltage	220 / 230 V AC
	Phase & Frequency	Single-phase 50 / 60 Hz
Weight	Gross weight	150 kg

• Dimensions

