

Slim!

PORTABLE X-RAY UNIT

AnyRad-5

- Features**
- Compact Size & Light Weight*
 - Superb & Stable Image Quality*
 - Flat Touch Panel & Digital Display*
 - Dose Area Product Function*
 - Operation by Mobile Apps*
 - 16 APR Data Selection*
 - Display Reverse Switch*
 - Built-in Dual Laser Pointer*
 - Easy access to DR Systems*







• Specifications

Model	AnyRad-5	
Generator	Type	High Frequency, Inverter
	Maximum Output	2.4 kW at 80 kV
	Frequency	40 kHz
	kV Range	40 - 100kV
	mA Range	16 - 40mA
	mAs Range	0.32 - 100mAs
	Indication	7-segment LED read out reversible
	APR	16 APR
X-ray Tube	Application	TOSHIBA (D-125S), Stationary
	Focal Spot	1.2 X 1.2mm
	Target Angle	16 Degree
	Anode Heat Storage	35 kJ[50kHU]
Collimator	Type	Manual
	Min. X-ray Field Size	5cm X 5cm at 100 cm SID
	Max. X-ray Field Size	35cm X 35cm at 66cm SID
	Lamp	Standard: Halogen, Option: LED
	Laser Pointer	Built in Dual Laser Pointer (Option)
Total Filtration	Inherent 0.8 mm Al, Added 2.0 mm Al	
Input Power rating	110/220 VAC(Free-Voltage), 50/60Hz, single phase	
Dimensions	485.6×244×192.5mm	
Weight	12.5Kg	
Standard Accessories	Hand Carry Aluminum Case	
	Power Cable (6M)	
	Hand Switch	
	Extra Halogen Lamp (1EA)	
	User Manual	

• Dimensions



Unit = mm