



# ANAESTHESIA

## LifeVet 8 Monitors

*Accurate and portable*



## Anaesthesia

# ACCURATE AND PORTABLE LifeVet 8M / 8C

---

*Our devices are designed with the needs of veterinary surgeons in mind. LifeVet 8 monitors are straightforward to use and provide reliable multi-parameter assessment during anaesthesia.*



**High  
Accuracy**

**Multiparameter  
Display**

### Rapid Documentation

The built in software allows paperless storage of data.

### Information Display

The high resolution screen allows clear visualisation of patients physiological parameters, providing continuous real time data.

### Enhanced patient safety

The optional ETCO<sub>2</sub> measurement feature via the side stream allows even greater levels of monitoring during anaesthesia.

### Compact Design

The machine's innovative design allows you to monitor patients in restricted spaces.

---

*“With the LifeVet 8 monitors I have peace of mind when my patients are under anaesthesia.”*


---

*Dependable*




## ANAESTHESIA LifeVet 8M/8C

- ▶ Pulse Oximetry with adjustable alarm and 360° visible alarm display
- ▶ Accurate respiratory rate, pulse rate and SpO<sub>2</sub> display
- ▶ 3-lead ECG for monitoring cardiac function
- ▶ Oesophageal or rectal temperature measurement
- ▶ Blood pressure monitoring (systolic/diastolic/mean arterial pressure (MAP))
- ▶ Capnometry (EtCO<sub>2</sub>) in the bypass measurement method (only LifeVet 8C)
- ▶ Mains or battery powered up to 4 hours



321900



321901

- ▶ High-resolution 8.4" LED touch screen (800 x 600 pixels)
- ▶ SpO<sub>2</sub> sensor with 2 sensor clips (small and large)
- ▶ Incl. universal set of alligator clips
- ▶ Incl. ECG cable, NiBP cable, extension cable and temperature probe
- ▶ Incl. Lithium ion battery
- ▶ Incl. blood pressure cuff for circumference: 6–11 cm
- ▶ Dimensions (in cm): W 27 x H 21 x D 11
- ▶ Weight (without battery): approx. 3.2 kg

### LifeVet 8M

**321900**

### LifeVet 8C

- ▶ Incl. CO<sub>2</sub> adapter, CO<sub>2</sub> line and water trap

**321901**

## Anaesthesia

# ACCESSORIES FOR LifeVet 8 Monitors

The following accessories support the use of LifeVet 8 monitors. All listed accessories are of the highest quality and have been tested for longevity.



### LINGUAL SENSOR FOR LIFEVET 8M / 8C

- ▶ For SpO<sub>2</sub> measurement
- ▶ Incl. two sensor clips (small and large) to accommodate different animal sizes

**321651**



### MOBILE STAND FOR LIFEVET 8M / 8C

- ▶ With basket and holding plate
- ▶ Secures the monitor with quick release

**321902**



### WALL BRACKET FOR LIFEVET 8M / 8C

**321903**



### BLOOD PRESSURE CUFF FOR LIFEVET 8M / 8C

- ▶ Circumference 3.1 – 5.7 cm

**E32181010**

- ▶ Circumference 4.3 – 8.0 cm

**321904**

- ▶ Circumference 5.8 – 10.9 cm

**321905**

- ▶ Circumference 7.1 – 13.1 cm

**E32181011**

- ▶ Circumference 8 – 15 cm

**E32181012**



### CO<sub>2</sub> WATER TRAP FOR LIFEVET 8C

**321906**

### CAPNOGRAPH SAMPLE LINE FOR LIFEVET 8C

- ▶ Length: 2.4 metres

- ▶ 10 per pack

**321378**

### CAPNOGRAPH CONNECTOR FOR LIFEVET 8C

- ▶ Inner/outer Ø: 15 mm

- ▶ Reusable

**321373**

### CO<sub>2</sub> AIRWAY ADAPTER FOR LIFEVET 8C

- ▶ Straight

- ▶ 10 per pack

**321907**